

PL-700 series Bench Top Meter

*pH · ORP · Conductivity · TDS · Salt · DO · Temperature
with Built-In Stirrer and USB Data Output*



Features

- Microprocessor based designed, and rugged housing with splash proof keyboard.
- The large backlight LCD displays pH or ORP or Conductivity or DO and Temperature simultaneously.
- Automatic Temperature Compensation(ATC), Manual Salinity(MSC) and Altitude(MAC) Compensation.
- Memory function stores and recalls up to 150 points. MAX/MIN and data Lock. Degree °C/°F are switchable.
- Combine special platform and electrode holder, or stirrer(Optional) for measuring convenient.
- The TDS value is determined by multiplying a conductivity reading by a known ratio in the range of 0.50 to 1.00.
- Output data for capturing on PC directly via USB.

